



PRESSALIT®

Product Item number	Pressalit Public Premium Set with waste bin, paper dispenser and soap dispenser, brushed steel PX30990
Product description	VVS: 773651315 NRF: 1438638 RSK: 8920373 EAN: 5708590392202 Wall mounting kit consisting of a one-hand operated 200 ml soap dispenser, paper dispenser for approx. 300 paper napkins and a 15 liter waste bin in brushed stainless steel. For wall mounting using glue or screws. The waste bin has an internal holder for standard waste sacks. The soap dispenser has a 200 ml container designed for regular liquid hand soap. The kit has a spring clip lock. Refilling of the paper dispenser is done from the top.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com