




PRESSALIT®

Product	Shower mattress for MSCT 1 and MSCT 4
Item number	R844717302
	VVS: 777542812 EAN: 5708590373904 HMI: 128437
Colour	Pewter grey
Product description	 Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 660 x 1660 mm. Incl. 1.6 m hose. Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 660 x 1660 mm. Incl. 1.6 m hose.
Materials	The mattress cover is welded polyester fabric with PVC-coating. It is UV-resistant and 100% impermeable to water. It contains neither CFC nor latex. The foam filling is flexible polyether-based polyurethane foam, approved in accordance with Oeko-tex 100, Class I.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.



PRESSALIT®

Product	Shower mattress for MSCT 1 and MSCT 4
Item number	R844717302
Cleaning	MSCT mattresses The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged. The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container. DISINFECTANTS: Sprit (husholdningssprit) Tvättssprit 70% Safe Zone Plus Selgiene Extreme Microsan Accel PREvention Wipes Lingettes The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.

pressalit.com

PRESSALIT®