




| | |
|----------------------------|---|
| Product | Mattress for MCT 2, MCT 3 and MCT 5 |
| Item number | R842717378 |
| | VVS: 777542718 EAN: 5708590373881 HMI: 128440 |
| Colour | Navy blue |
| Product description |  <p>Mattress for MCT 2, MCT 3 and MCT 5 change trolleys, 695 x 1680 mm. Can also be used with MSCT 1</p> |
| Materials | The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I. |
| Warranty | Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date. |
| Cleaning | <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p>DISINFECTANTS: Sprit (husholdningsprit) Tvättsprit 70% Safe Zone Plus Selgiene Extreme Accel PREvention Wipes Lingettes OxiGenesis from ENVIROX</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p> |



PRESSALIT®

| | |
|--------------------|-------------------------------------|
| Product | Mattress for MCT 2, MCT 3 and MCT 5 |
| Item number | R842717378 |

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®