



PRESSALIT®

Product Item number	Pressalit Choice Spare paper holder, polished steel Q10260
Product description	VVS:776591240 NRF:6190933 RSK:8744046 EAN:5708590352831  Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®