



|                            |  |
|----------------------------|--|
| <b>Product</b>             | Pressalit Choice Toilet brush for wall, brushed brass  |
| <b>Item number</b>         | Q70850   |
| <b>Product description</b> | <p>VVS:778497205<br/>NRF:7001083<br/>RSK:8714275<br/>EAN:5708590372242<br/>DB: 2173273</p> <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p> <p>Elegant toilet brush for wall mounting. Opens up extra space on the floor.</p> |
| <b>Materials</b>           | Stainless steel AISI 304, mounting knot in brass   |
| <b>Warranty</b>            | A 2-year warranty is provided for material and manufacturing defects.  |
| <b>Cleaning</b>            | <p>This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>  |