




PRESSALIT®

Product Item number	Pressalit Public Classic Waste bin 60 l, matte black PY40920
Product description	VVS:777744631 NRF:1438612 RSK:8920402 EAN:5708590392417  Stainless steel waste bin for floor mounting or concealed wall mounting. Front panel with a stainless steel plate with rounded corners. Equipped with removable, washable polypropylene fabric bag. Capacity approx. 60 liters. Comes with mounting kit. 2.5 mm and 4 mm thick grating and mounting holes placed in the second row for concealed wall mounting.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com