




PRESSALIT®

Product	Pressalit Trend incl. Pressalit Sway D2 994
Item number	CP994011CF
	VVS:613326500 NRF:1462969 RSK:7822116 EAN:5708590396903
Colour	White
Product description	The Pressalit Trend Combi Pack with Sway D2 seat offers a perfect union of minimalist elegance and advanced functionality. The rimless, wall-hung toilet features the hygienic and efficient Circular Flush system for a quiet, splash-free flush. The award-winning Sway D2 seat is slim and stylish, complementing the bowl's modern design while providing excellent comfort. Soft Close ensures a gentle, noiseless close, and Lift-Off makes cleaning simple. An antibacterial additive helps reduce bacterial growth for improved hygiene. This combination is ideal for contemporary bathrooms where design and practicality are equally important, creating a refined, easy-to-maintain solution for everyday use.
Product capacity	Load - ring seat 240 kg





PRESSALIT®

Product	Pressalit Trend incl. Pressalit Sway D2 994
Item number	CP994011CF
Materials	BOWL: Porcelain TOILET SEAT: The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises approx. 70% urea formaldehyde resin, 30% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. Free formaldehyde is released during the production process. Pressalit toilet seats have been measured after curing ($< 0.01 \text{ mg/m}^3$) to emit significantly below the permitted EU limit of 0.08 mg/m^3. ANTIBACTERIAL: This toilet seat features an antibacterial additive with biocidal silver phosphate glass, which helps prevent bacterial growth. BUFFERS: EVA (copolymer made of ethylene and vinyl acetate). DAMPER KIT: Hydraulic damper with damper housing in stainless steel or thermoplastic, mounted with plastic parts in thermoplastic. MOUNT/HINGES: Stainless steel and thermoplastic.
Warranty	Pressalit provides a 10-year warranty on toilet seats and a 5-year warranty on toilet bowls against presentation of a valid receipt. In the event of a complaint, please contact the retailer or store where the product was purchased. The warranty covers defective materials and/or manufacturing defects and/or functional limitations that may occur following proper use and cleaning of the product. Pressalit reserves the right to examine the product at our factory. In the case of a valid complaint, Pressalit will either: - supply spare parts, or - replace the product with an identical or a similar model. If the color of the toilet seat is no longer available, Pressalit reserves the right to supply a similar model in white. Upon replacement or delivery of spare parts, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit is entitled instead to refund the purchase price. Subject to changes in the product range, technical changes and printing errors.



PRESSALIT®

Product	Pressalit Trend incl. Pressalit Sway D2 994
Item number	CP994011CF

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

23. september 2026