




PRESSALIT®

Product Item number	Pressalit Public Premium Toilet roll holder for 1 large roll of toilet paper, matte black PX40225
Product description	VVS: 776130131 RSK: 8920376 EAN: 5708590392714  <p>Toilet paper holder in large steel housing with space for rolls of max. Ø 280 x 100 mm. For wall mounting. The holder is lockable and includes mounting kit.</p>
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®