




PRESSALIT®

Product	Wall-hung toilet 700 mm
Item number	R2099000
Colour	White
Product description	<p>VVS:612928000 NRF:6191026 RSK:7822154 EAN:5708590345284 HMI: 105700</p> <p>Toilet, wall-hung, distance from wall to the front: 700 mm. Excl. mounting kit.</p> <p>Toilet, wall-hung, distance from wall to the front: 700 mm. Excl. mounting kit.</p>
Product capacity	Max load: 400 kg.
Materials	Sanitary porcelain
Ergonomy	700 mm is recommended for wheelchair users as it gives more space for a carer. Wall-hung toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a mobile hoist, and the carer's feet.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®