



PRESSALIT®

Product Item number	Pressalit Public Universal Soap dispenser with sensor, 850 ml, brushed steel PZ30715
Product description	VVS: 774566215 NRF: 1438658 RSK: 8920590 EAN: 5708590396347 Soap dispenser with stainless steel sensor for wall mounting. Drip-free soap pump with 850 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®