




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Classic Toilet brush for wall mounting, w/ bowl, brushed steel PY30810
<b>Product description</b>	<p>VVS: 778499115 NRF: 1438655 RSK: 8920393 EAN: 5708590392356</p>  <p>Toilet brush in solid stainless steel for wall mounting. Concealed anti-theft protection with 2-point rosette. Incl. mounting kit.</p>
<b>Materials</b>	European stainless steel (AISI 303)
<b>Warranty</b> <b>Cleaning</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

[pressalit.com](http://pressalit.com)

## PRESSALIT®