



PRESSALIT®

Splashguard, transparent mat, for 50 mm raised toilet seats

Item number

T9516
VVS: 616810516
NRF: 6190363
EAN: 5708590183909

Description

Splashguard, transparent mat, for 50 mm raised toilet seats..

Splashguard, transparent mat, for 50 mm raised toilet seats.

**Materials**

Polypropylene (PP)

Product capacity**Weight**

0,5 kg

Gross weight

0,5 kg

Weight per crate

0 kg

Qty per crate

0

Qty per pallet

0

pressalit.com
