



PRESSALIT®

Splash guard, matt transparent, for 100 mm raised toilet seats

Item number

T9515
VVS: 616810515
NRF: 6190362
EAN: 5708590070551

Description

Splash guard, matt transparent, for 100 mm raised toilet seats



Materials

Polypropylene (PP)

Product capacity**Weight**

0,5 kg

Gross weight

0,5 kg

Weight per crate

0 kg

Qty per crate

0

Qty per pallet

0

The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

pressalit.com
