




PRESSALIT®

Product	Safety lock for CT 4000 manually foldable change table
Item number	R8486
	VVS: 777545983 EAN: 5708590371191
Product description	 The safety lock holds the change table securely folded against the wall. Safety lock, 330 mm, wall-mounted, for CT 4000 manually foldable change table

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®