



PRESSALIT®

Product Item number	Pressalit Public Classic Soap dispenser with sensor, 850 ml, polished steel PY10715
Product description	VVS: 774565240 NRF: 1438641 RSK: 8920382 EAN: 5708590392264 Soap dispenser with stainless steel sensor for wall mounting. Drip-free soap pump with 850 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com