




PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, brushed steel PY30725
Product description	VVS: 774565615 NRF: 1438614 RSK: 8920390 EAN: 5708590392745  Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning WEEE	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Please enquire of your local Pressalit subsidiary or distributor (in countries where there is no subsidiary) that will arrange for recycling of electronic waste in accordance with the WEEE directive. Subject to changes in the product range, technical changes and printing errors.



PRESSALIT®

Product	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, brushed steel
Item number	PY30725

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

24. september 2026