



PRESSALIT®

Product	Wireless remote for SCT 1000/2000/2100/3000 and CT 4000
Item number	R8493
	EAN:5708590378350
Product description	Can be purchased as additional accessory or as a spare part. Wireless remote for SCT 1000/2000/2100/3000 and CT 4000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®