



# PRESSALIT®

---

## Raised buffer with screw 50 mm

---


**Item number**

T9803

EAN:5708590183947

**Description**

Raised buffer with screw 50 mm



---

**Product capacity**

**Weight**

0,46 kg

**Gross weight**

0,46 kg

**Weight per crate**

0 kg

**Qty per crate**

0

**Qty per pallet**

0

The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

---

[pressalit.com](http://pressalit.com)

---