



<b>Product</b> <b>Item number</b>	Pressalit Choice Towel hook, curve, 2 pcs., brushed brass Q70160
<b>Product description</b>	VVS:775295305 NRF:7001574 RSK:8920049 EAN:5708590381862 DB: 2425188  Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.  Use the towel hook in the bathroom and in the kitchen. That way, you can let the style be consistent throughout the house.
<b>Materials</b>	Stainless steel AISI 304, mounting knot in brass
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.