



# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Classic Soap dispenser with sensor, 850 ml, brushed steel PY30715
<b>Product description</b>	VVS:774565215 NRF:1438639 RSK:8920389 EAN:5708590392738  Soap dispenser with stainless steel sensor for wall mounting. Drip-free soap pump with 850 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)