



PRESSALIT®

Product Item number	Pressalit Style Towel rack, double, 31.89", chrome PR20560
Product description	VVS:775712404 NRF:7001529 RSK:8755168 EAN:5708590375205 DB: 2266849 Let elegance and your personal style shine through with Pressalit Style. The long, elegant, graphic shapes will lend an exciting touch to your bathroom. Not only will your bathroom get an aesthetic makeover - it will also reflect the fact that you have made an active decision to welcome extraordinary Danish design into the décor of your home. The products are easy to fit - with either included screws or glue, which is purchased separately.
Materials	Stainless steel AISI 304 and Zamak
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®