



# PRESSALIT®

## Pressalit Trend incl. Pressalit Comfort D2 1050

Wall hung toilet with circular flush, incl. toilet seat  
Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel

### Item number




CP1050011CF  
VVS: 613326400  
NRF: 1462954  
RSK: 7822115  
EAN: 5708590396897

The Pressalit Trend Combi Pack with Comfort D2 seat combines a stylish wall-hung toilet bowl with a perfectly matched seat for exceptional comfort and stability. The rimless bowl with innovative Circular Flush ensures powerful, splash-free and quiet flushing, while the wrap-over Comfort D2 seat blends seamlessly with the design. Pressalit SteadyFix® mounting provides outstanding stability, and Soft Close plus Lift-Off functions make cleaning easy and use effortless. The antibacterial additive helps inhibit bacterial growth. This durable, high-quality combination delivers elegance, functionality, and hygiene in perfect harmony. A complete bathroom solution that is both practical and beautifully designed.



### Colour

White





---

**Pressalit Trend incl. Pressalit Comfort D2 1050**

**Wall hung toilet with circular flush, incl. toilet seat<br><br>Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel**

---

**Materials****TOILET SEAT:**

The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises 67% urea formaldehyde resin, 28% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. The raw material from the supplier contains up to 0.03% free formaldehyde and even less following curing.

**ANTIBACTERIAL:**

This toilet seat features an antibacterial additive with biocidal silver chloride, which helps prevent bacterial growth.

**BUFFERS:**

EVA (copolymer made of ethylene and vinyl acetate).

**DAMPER KIT (only relevant for toilet seats with soft close):**

Hydraulic damper with stainless steel damper casing and thermoplastic parts.

**MOUNT/HINGES:**

Stainless steel model W. no. 1.4301/M and thermoplastic.

**Product capacity**

Load - ring seat 240 kg

**Weight**

31,034 kg

**Carton size**

540x410x510 mm

**Gross weight**

33,703 kg

**Weight per crate**

33,703 kg

**Qty per crate**

1

**Crates per pallet**

5

**Qty per pallet**

5

**Warranty**

Pressalit provides a 10-year warranty on toilet seats and a 5-year warranty on toilet bowls against presentation of a valid receipt. In the event of a complaint, please contact the retailer or store where the product was purchased. The warranty covers defective materials and/or manufacturing defects and/or functional limitations that may occur following proper use and cleaning of the product. Pressalit reserves the right to examine the product at our factory. In the case of a valid complaint, Pressalit will either:<br>- supply spare parts, or<br>- replace the product with an identical or a similar model.<br>If the color of the toilet seat is no longer available, Pressalit reserves the right to supply a similar model in white. Upon replacement or delivery of spare parts, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit is entitled instead to refund the purchase price.



# PRESSALIT®

---

**Pressalit Trend incl. Pressalit Comfort D2 1050**

**Wall hung toilet with circular flush, incl. toilet seat<br><br>Toilet  
seat with soft close, lift-off and antibacterial effect, incl. hinge in  
stainless steel**

---

Subject to changes in the product range, technical changes and printing errors.

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®