



PRESSALIT®

Product Item number	Mounting kit V8666 (8 pcs), for aerated concrete V8666
Product description	EAN:5708590351285 Mounting kit, 8 bolts M5 x 50, plugs FTP K 6, for aerated concrete. Used for: PLUS grab bars: RT018, RT024, RT036, RT042 PLUS modular shelves: RT490, RT491 Mounting kit, 8 bolts M5 x 50, plugs FTP K 6, for aerated concrete
Materials	Stainless steel and plastics

pressalit.com