



PRESSALIT®

Pressalit Trend incl. Pressalit Sway D2 994
Wall hung toilet with circular flush, incl. toilet seat
seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel

Item number


CP994011CF
VVS: 613326500
NRF: 1462969
RSK: 7822116
EAN: 5708590396903

The Pressalit Trend Combi Pack with Sway D2 seat offers a perfect union of minimalist elegance and advanced functionality. The rimless, wall-hung toilet features the hygienic and efficient Circular Flush system for a quiet, splash-free flush. The award-winning Sway D2 seat is slim and stylish, complementing the bowl's modern design while providing excellent comfort. Soft Close ensures a gentle, noiseless close, and Lift-Off makes cleaning simple. An antibacterial additive helps reduce bacterial growth for improved hygiene. This combination is ideal for contemporary bathrooms where design and practicality are equally important, creating a refined, easy-to-maintain solution for everyday use.




Colour

White



*20 mm / ** 40 mm





Pressalit Trend incl. Pressalit Sway D2 994
Wall hung toilet with circular flush, incl. toilet seat
seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel

Materials

BOWL:
Porcelain

TOILET SEAT:

The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises approx. 70% urea formaldehyde resin, 30% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. Free formaldehyde is released during the production process. Pressalit toilet seats have been measured after curing (< 0.01 mg/m³) to emit significantly below the permitted EU limit of 0.08 mg/m³.

ANTIBACTERIAL:

This toilet seat features an antibacterial additive with biocidal silver phosphate glass, which helps prevent bacterial growth.

BUFFERS:

EVA (copolymer made of ethylene and vinyl acetate).

DAMPER KIT:

Hydraulic damper with damper housing in stainless steel or thermoplastic, mounted with plastic parts in thermoplastic.

MOUNT/HINGES:

Stainless steel and thermoplastic.

Product capacity

Load - ring seat 240 kg

Weight

32,035 kg

Carton size

540x410x510 mm

Gross weight

34,704 kg

Weight per crate

34,704 kg

Qty per crate

1

Crates per pallet

5

Qty per pallet

5



PRESSALIT®

Pressalit Trend incl. Pressalit Sway D2 994
Wall hung toilet with circular flush, incl. toilet seat
seat with soft close, lift-off and antibacterial effect, incl. hinge in
stainless steel

Warranty

Pressalit provides a 10-year warranty on toilet seats and a 5-year warranty on toilet bowls against presentation of a valid receipt. In the event of a complaint, please contact the retailer or store where the product was purchased. The warranty covers defective materials and/or manufacturing defects and/or functional limitations that may occur following proper use and cleaning of the product. Pressalit reserves the right to examine the product at our factory. In the case of a valid complaint, Pressalit will either:
- supply spare parts, or
- replace the product with an identical or a similar model.
If the color of the toilet seat is no longer available, Pressalit reserves the right to supply a similar model in white. Upon replacement or delivery of spare parts, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit is entitled instead to refund the purchase price.

Subject to changes in the product range, technical changes and printing errors.

pressalit.com
