



PRESSALIT®

Product Item number	Pressalit Public Classic Soap dispenser, 1100 ml, matt black PY40710
Product description	VVS: 774565131 NRF: 1438643 RSK: 8920397 EAN: 5708590392394 Stainless steel soap dispenser for wall mounting. One-handed, drip-free soap pump with 1100 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com