



<b>Product</b> <b>Item number</b>	Pressalit Public Classic Soap dispenser, 1100 ml, polished steel PY10710
<b>Product description</b>	VVS:774565140 NRF:1438644 RSK:8920381 EAN:5708590392257  Stainless steel soap dispenser for wall mounting. One-handed, drip-free soap pump with 1100 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)