




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Premium Toilet paper holder for 1 roll of toilet paper + 1 spare roll, brushed steel PX30210
<b>Product description</b>	<p>VVS: 776130215 NRF: 1438646 RSK: 8920365 EAN: 5708590392660</p>  <p>Toilet paper holder in steel housing with space for spare roll. For wall mounting. The holder is lockable and includes mounting kit.</p>
<b>Materials</b>	European stainless steel (AISI 304)
<b>Warranty</b> <b>Cleaning</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

[pressalit.com](http://pressalit.com)

# PRESSALIT®