




PRESSALIT®

Universal hinge

Item number

DK4999
VVS: 616810820
EAN: 5708590370057



Product capacity

Weight

0,146 kg

Gross weight

0,146 kg

Weight per crate

0 kg

Qty per crate

0

Qty per pallet

0

The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®