



PRESSALIT®

Pressalit Trend incl. Pressalit Comfort D2 1050

Wall hung toilet with circular flush, incl. toilet seat
Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel

Item number




CP1050011CF
VVS: 613326400
NRF: 1462954
RSK: 7822115
EAN: 5708590396897

The Pressalit Trend Combi Pack with Comfort D2 seat combines a stylish wall-hung toilet bowl with a perfectly matched seat for exceptional comfort and stability. The rimless bowl with innovative Circular Flush ensures powerful, splash-free and quiet flushing, while the wrap-over Comfort D2 seat blends seamlessly with the design. Pressalit SteadyFix® mounting provides outstanding stability, and Soft Close plus Lift-Off functions make cleaning easy and use effortless. The antibacterial additive helps inhibit bacterial growth. This durable, high-quality combination delivers elegance, functionality, and hygiene in perfect harmony. A complete bathroom solution that is both practical and beautifully designed.



Colour

White





PRESSALIT®

Pressalit Trend incl. Pressalit Comfort D2 1050

Wall hung toilet with circular flush, incl. toilet seat
Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel

Materials

BOWL:
Porcelain

TOILET SEAT:

The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises approx. 70% urea formaldehyde resin, 30% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. Free formaldehyde is released during the production process. Pressalit toilet seats have been measured after curing ($< 0.01 \text{ mg/m}^3$) to emit significantly below the permitted EU limit of 0.08 mg/m^3 .

ANTIBACTERIAL:

This toilet seat features an antibacterial additive with biocidal silver phosphate glass, which helps prevent bacterial growth.

BUFFERS:

EVA (copolymer made of ethylene and vinyl acetate).

DAMPER KIT:

Hydraulic damper with damper housing in stainless steel or thermoplastic, mounted with plastic parts in thermoplastic.

MOUNT/HINGES:

Stainless steel and thermoplastic.

Product capacity

Load - ring seat 240 kg

Weight

31,034 kg

Carton size

540x410x510 mm

Gross weight

33,703 kg

Weight per crate

33,703 kg

Qty per crate

1

Crates per pallet

5

Qty per pallet

5



PRESSALIT®

Pressalit Trend incl. Pressalit Comfort D2 1050

**Wall hung toilet with circular flush, incl. toilet seat

Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel**

Warranty

Pressalit provides a 10-year warranty on toilet seats and a 5-year warranty on toilet bowls against presentation of a valid receipt. In the event of a complaint, please contact the retailer or store where the product was purchased. The warranty covers defective materials and/or manufacturing defects and/or functional limitations that may occur following proper use and cleaning of the product. Pressalit reserves the right to examine the product at our factory. In the case of a valid complaint, Pressalit will either:
- supply spare parts, or
- replace the product with an identical or a similar model.
If the color of the toilet seat is no longer available, Pressalit reserves the right to supply a similar model in white. Upon replacement or delivery of spare parts, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit is entitled instead to refund the purchase price.

Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®