



Product	Toilet paper holder for TMA fold down grab bars
Item number	R1175000
Color	White
Product description	<p>VVS:615723950 NRF:6191655 RSK:8920763 EAN:5708590211725 HMI: 35271 NBS:N13:436 DB: 2076625</p> <p>Paper holder for TMA arm supports. Fixed, can be mounted on the left or right handle.</p> <p>Paper holder for TMA arm supports. Fixed, can be mounted on the left or right handle.</p>
Adjustability	Fixed mounting
Materials	Powder-coated stainless steel.
Ergonomy	Can be mounted at the right or the left handle.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.

