




PRESSALIT®

Product Item number	Mounting kit V8654 (4 pcs), for aerated concrete V8654
Product description	EAN:5708590351407  Mounting kit, 4 countersunk bolts M8 x 35, plugs FTP M 8, for aerated concrete. Used for: PLUS support arms: R360070, R360085 PLUS folding seat: R5510, R5504 PLUS shower seats: R7304, R7404 Mounting kit, 4 countersunk bolts M8 x 35, plugs FTP M 8, for aerated concrete
Materials	Stainless steel and plastic

pressalit.com

PRESSALIT®