




PRESSALIT®

Product Item number	Freestanding column 37.2", for height adjustable PLUS fold down grab bar R9431
Product description	VVS:615856990 NRF:6190668 EAN:5708590320434 NBS:N13:436  Free-standing column, gray, for height adjustable PLUS fold down grab bar Free-standing column, gray, for height adjustable PLUS support arm
Adjustment height	9,46 inch
Materials	Core material: Steel Surface material: Rilsan
Features	Used with freestanding toilets and when it is not possible to wallmount a grab bar. For height adjustable fold down grab bars.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.

pressalit.com

PRESSALIT®