

DECLARATION OF PERFORMANCE

No. CombiPack-EN997

1. Unique identification code of the product-type:

Product type: Combi Pack

For product number please see the attached list.

2. Intended use/es:

Personal hygiene.

3. Manufacturer:

Pressalit A/S

Pressalit 1

8660 Ry

4. System/s of AVCP:

System 4

5. Harmonised standard:

EN 997:2018

6. Declared performance/s:

Essential characteristics	Performance	Harmonised technical specification
Capacity of flushing water (CF)	Pass	EN 997:2018
Backflow prevention (of foul air): (class 1) (BP)	Pass	
Cleanability (CA)	Pass	
Load resistance (LR)	Pass	
Water/leak tightness (WL)	NPD*	
Valve reliability (VR)	NPD*	
Durability (DA)	Pass	

* No performance determined

The performance of the product identified above is in conformity with the set of declared performance/s. This declaration of performance is issued, in accordance with Regulation (EU) No 305/2011, under the sole responsibility of the manufacturer identified above.

Signed for and on behalf of the manufacturer by:

Date of issue: 26.05.2026

At: Ry



Pressalit A/S, Mette Elstrøm Kristoffersen, Regulatory Specialist

Attachment to Declaration of Performance

No. CombiPack-EN997

Product no.	WC Pan	Name	Designation code
CP1050011CF	A10279	Pressalit Trend wall hung toilet with circular flush incl. toilet seat Pressalit Comfort D2 with soft close and lift-off, antibacterial, white.	EN 997 - CL 1 - 6 A
CP1008011CF	A10279	Pressalit Trend wall hung toilet with circular flush incl. toilet seat Projecta D Solid Pro with soft close, antibacterial, white.	EN 997 - CL 1 - 6 A
CP1206011CF	A10279	Pressalit Trend wall hung toilet with circular flush, incl. toilet seat Pressalit Style D2 with soft close and lift-off, antibacterial, white	EN 997 - CL 1 - 6 A
220366	A10279	Pressalit Trend wallhung toilet with circular flush, white	EN 997 - CL 1 - 6 A
CP1206111CF	A10279	Pressalit Trend wall hung toilet with circular flush, incl. toilet seat Pressalit Style D2 with soft close and lift-off, antibacterial, black	EN 997 - CL 1 - 6 A
CP1206230CF	A10279	Pressalit Trend wall hung toilet with circular flush, incl. toilet seat Pressalit Style D2 with soft close and lift-off, matt white	EN 997 - CL 1 - 6 A
CP1206231CF	A10279	Pressalit Trend wall hung toilet with circular flush, incl. toilet seat Pressalit Style D2 with soft close and lift-off, matt black	EN 997 - CL 1 - 6 A
CP1206367CF	A10279	Pressalit Trend wall hung toilet with circular flush, incl. toilet seat Pressalit Style D2 with soft close and lift-off, matt dusty green	EN 997 - CL 1 - 6 A
CP1206085CF	A10279	Pressalit Trend wall hung toilet with circular flush, incl. toilet seat Pressalit Style D2 with soft close and lift-off, matt taupe	EN 997 - CL 1 - 6 A
CP994011CF	A10279	Pressalit Trend wall hung toilet with circular flush, white incl. toilet seat Pressalit Sway D2, with soft close and lift-off, antibacterial, white	EN 997 - CL 1 - 6 A
CP994111CF	A10279	Pressalit Trend wall hung toilet with circular flush, white incl. toilet seat Pressalit Sway D2, with soft close and lift-off, antibacteriel, black	EN 997 - CL 1 - 6 A
CP994159CF	A10279	Pressalit Trend wall hung toilet with circular flush, white incl. toilet seat Pressalit Sway D2, with soft close and lift-off, matt royal blue	EN 997 - CL 1 - 6 A
CP994230CF	A10279	Pressalit Trend wall hung toilet with circular flush, white incl. toilet seat Pressalit Sway D2, with soft close and lift-off, matt white	EN 997 - CL 1 - 6 A
CP994231CF	A10279	Pressalit Trend wall hung toilet with circular flush, white incl. toilet seat Pressalit Sway D2, with soft close and lift-off, matt black	EN 997 - CL 1 - 6 A