



PRESSALIT®

Product	Mounting kit V8666 (8 pcs), for aerated concrete
Item number	V8666
	EAN:5708590351285
Product description	Mounting kit, 8 bolts M5 x 50, plugs FTP K 6, for aerated concrete. Used for: PLUS grab bars: RT018, RT024, RT036, RT042 PLUS modular shelves: RT490, RT491 Mounting kit, 8 bolts M5 x 50, plugs FTP K 6, for aerated concrete
Materials	Stainless steel and plastics

pressalit.com

PRESSALIT®