




PRESSALIT®

Product Item number	Pressalit Public Classic Soap dispenser with sensor, 850 ml, matt black PY40715
Product description	VVS: 774565231 RSK: 8920398 EAN: 5708590392752  <p>Soap dispenser with stainless steel sensor for wall mounting. Drip-free soap pump with 850 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.</p>
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®