




PRESSALIT®

Product	Shower mattress for MSCT 1 and MSCT 4
Item number	R844717302
	VVS: 777542812 EAN: 5708590373904 HMI: 128437
Color	Pewter gray
Product description	 <p>Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 26" x 65.4". Incl. 63" hose.</p> <p>Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 26" x 65.4". Incl. 63" hose.</p>
Materials	The mattress cover is welded polyester fabric with PVC-coating. It is UV-resistant and 100% impermeable to water. It contains neither CFC nor latex. The foam filling is flexible polyether-based polyurethane foam, approved in accordance with Oeko-tex 100, Class I.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.
Cleaning	<p>This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p>DISINFECTANTS: Sprit (husholdningsprit) Tvättsprit 70% Safe Zone Plus Selgiene Extreme Microsan Accel PREvention Wipes Lingettes</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>



PRESSALIT®

Product	Shower mattress for MSCT 1 and MSCT 4
Item number	R844717302

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®