



# United Kingdom Conformity Assessment

**United Kingdom RP:**  
 Pressalit Ltd.  
 100, Longwater Avenue  
 Green Park  
 Reading RG2 6GP  
 United Kingdom  
 T: 01180923 2696  
 uk@pressalit.com

This product is produced in compliance with all relevant requirements of:

- Regulation (EU) 2017/745 on medical devices (Risk class I, rule 1)
- Council Directive 2006/42/EC on machinery
- Council Directive 2011/65/EC on RoHS

For the UK:

- Medical Device Regulation (MDR) 2002 on medical devices
- Supply of Machinery (Safety) Regulations 2008
- Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

Medical devices - Application of risk management and quality management in design and production

ISO 14971:2019  
 ISO 13485:2016

The products are manufactured in accordance with the following standards:

IEC 60601-1:2005  
 ISO 17966:2016

Ry, 23. juni 2026

**Mette Elstrøm Kristoffersen, QA/RA Specialist, Pressalit A/S**  
 Responsible

**Manufacturer:** Pressalit A/S - Pressalitvej 1 DK-8580 Ry - T: +45 8788 8788 - pressalit@pressalit.com

## SELECT TL5



**R9221000**

**EAN (GTIN):** 5708590390963

**Nomenclature (GMDN):** 63055

**Nomenclature (EMDN):** Y091299

### Name

EN SELECT TL5 toilet lifter with tilt function, electrical, for wall outlet

### Intended purpose

EN Toilet lifters ensure optimum ergonomics at a toilet, for both seated and standing users with disabilities. Expected service life is 10 years, provided that maintenance is carried out as described in the user manual.