




Product Item number	TMA 2 arm support for mounting on toilet, wider fit R1191112000
Colour Product description	VVS:615724100 NRF:6190325 EAN:5708590376240 HMI: 130727 DB: 2510918 FDA:ILS/D135866 White  <p>TMA arm supports comes in various widths, and models with or without toilet seat.</p> <p>TMA 2 arm supports for mounting on toilet, incl. BL1999 hinge. Dark grey handles. BL1 hinge, is only recommended with toilet seats R26, R27, R28, R29, R37 or R38.</p>
Adjustability	Fixed mounting
Materials Ergonomy Features	Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i> The arm supports provide firm support and the moulded handles are comfortable for the hands and arms. The arm supports can be folded up and down individually.
Warranty Cleaning	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged. The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.



PRESSALIT®

Product	TMA 2 arm support for mounting on toilet, wider fit
Item number	R1191112000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®