




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Choice Toilet brush for wall with glass bowl, polished steel Q10800
<b>Product description</b>	<p>VVS:778497140 NRF:6190942 RSK:8714259 EAN:5708590336947</p>  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p> <p>Elegant toilet brush for wall mounting. Opens up extra space on the floor.</p>
<b>Materials</b>	Stainless steel AISI 304, glass, mounting knot in brass
<b>Warranty</b>	A 2-year warranty is provided for material and manufacturing defects.
<b>Cleaning</b>	<p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

[pressalit.com](http://pressalit.com)