




PRESSALIT®

| | |
|--------------------------------------|---|
| Product Item number | Freestanding column 945 mm, for height adjustable PLUS support arm R9431 |
| Product description | <p> VVS:615856990 NRF:6190668 EAN:5708590320434 NBS:N13:436 </p>  <p> Free-standing counn, grey, for height adjustable PLUS support arm Free-standing counn, grey, for height adjustable PLUS support arm </p> |
| Adjustment height | 240 mm |
| Materials | Basic material: Steel Surface material: Rilsan |
| Features | Used with freestanding toilets and where it is not possible to wallmount a support arm. For adjustable height support arm. |
| Warranty | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date. |



PRESSALIT®

| | |
|--------------------|--|
| Product | Freestanding column 945 mm, for height adjustable PLUS support arm |
| Item number | R9431 |

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®