




Product	Changing table mattress 1350 mm
Item number	R8705000
Colour	White
Product description	<p>VVS:777630830 RSK:8565138 EAN:5708590366111 HMI: 125814</p>  <p>Mattresses for daycare changing tables are available in several colors and in lengths that are adapted to the different changing table lengths.</p> <p>Accessories for changing tables. Mattress 760 mm x 1350 mm x 30 mm, for 1800 mm changing table with sanitary appliances. Phthalate-free</p>
Materials	Cover: Techmadin, phthalate-free Padding: Foam
Cleaning	<p>Upholstery (Cover): Techmadin, phthalate-free Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and wash regularly at 60-90°, without fabric softener. Shrinkage after washing may occur. Do not turn inside out before washing. Can be tumble-dried at max. 60°. Can be disinfected with antiseptic wipes. Alcohol % < 70. Do not use disinfectants intended for hard surfaces.</p> <p>Padding: (foam) The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.</p>



PRESSALIT®

Product	Changing table mattress 1350 mm
Item number	R8705000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®