




PRESSALIT®

<b>Product</b>	Changing table mattress 1400 mm
<b>Item number</b>	R8702000
<b>Colour</b>	White
<b>Product description</b>	<p>VVS:777630730 RSK:8565127 EAN:5708590366067 HMI: 125781</p>  <p>Mattresses for daycare changing tables are available in several colors and in lengths that are adapted to the different changing table lengths.</p> <p>Accessories for changing tables Mattress 760 mm x 1400 mm x 30 mm, for 1400 mm changing table Phthalate-free</p>
<b>Materials</b>	Cover: Techmadin, phthalate-free Padding: Foam
<b>Cleaning</b>	<p><b>Upholstery (Cover):</b> Techmadin, phthalate-free Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and wash regularly at 60-90°, without fabric softener. Shrinkage after washing may occur. Do not turn inside out before washing. Can be tumble-dried at max. 60°. Can be disinfected with antiseptic wipes. Alcohol % &lt; 70. Do not use disinfectants intended for hard surfaces.</p> <p><b>Padding: (foam)</b> The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.</p>



# PRESSALIT®

<b>Product</b>	Changing table mattress 1400 mm
<b>Item number</b>	R8702000

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

PRESSALIT®