



PRESSALIT®

Product Item number	Pressalit Style Toilet wastebasket, brushed brass/black PR70950
Product description	VVS:777742105 NRF:7001571 RSK:8900162 EAN:5708590377575 DB: 2266886 Let elegance and your personal style shine through with Pressalit Style. The long, elegant, graphic shapes will lend an exciting touch to your bathroom. Not only will your bathroom get an aesthetic makeover - it will also reflect the fact that you have made an active decision to welcome extraordinary Danish design into the décor of your home. The products are easy to fit - with either screws or glue, which is purchased separately. Complete your bathroom décor with elegant details. Removable 5.0-liter inner bucket. The toilet waste basket opens with a light touch. You can stand it on the floor or hang it on the wall if you wish to make best use of your floor space and facilitate cleaning. The discreet rubber ring on the bottom ensures greater stability if you choose to place the toilet waste basket on the floor.
Materials	Stainless steel AISI 304, ABS, PP, TPE and PA
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.



PRESSALIT®

Product	Pressalit Style Toilet wastebasket, brushed brass/black
Item number	PR70950

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

24. september 2026