




PRESSALIT®

Product Item number	Pressalit Public Premium Toilet paper holder with space for 5 rolls of toilet paper, brushed steel PX30220	
Product description	VVS: 776130515 NRF: 1438652 RSK: 8920366 EAN: 5708590392141	
	Stainless steel toilet paper holder with space for storing 5 rolls. For wall mounting. Designed for one standard toilet roll and four hidden rolls that automatically drop down. Including mounting kit. Inspection opening to indicate filling level. Refill through opening in the top.	
Materials	European stainless steel (AISI 304)	
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.	

pressalit.com