




PRESSALIT®

Product Item number	TMA 1 arm support for mounting on toilet R1171000
Colour	VVS:615723100 NRF:6191651 EAN:5708590211718 HMI: 42230 NBS:N13:436 DB: 2076624 FDA:ILS/D135866 White
Product description	 TMA arm supports comes in various widths, and models with or without toilet seat. TMA 1 arm supports for mounting on toilet, incl. BL1999 hinge. Dark grey handles. BL1 hinge, is only recommended with toilet seats R26, R27, R28, R29, R37 or R38.
Adjustability	Fixed mounting
Materials	Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i>
Ergonomy	The folding arms give a good support and the moulded handles are comfortable for the hands and arms.
Features	The arm supports can be folded individually.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
Cleaning	The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged. The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.



PRESSALIT®

Product	TMA 1 arm support for mounting on toilet
Item number	R1171000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

8. juli 2026