




|                                                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|--------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b><br><b>Item number</b>                   | VALUE IV folding shower seat with supporting leg, fixed height<br>R1622000112                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Colour</b>                                          | VVS:777525412<br>NRF:6314098<br>RSK:8565070<br>EAN:5708590377728<br>HMI: 124486<br>DB: 2510921<br>FDA:INN/D135865<br>Anthracite grey                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Product description</b>                             |  <p>VALUE IV folding shower seat and backrest create a simple solution for users who need to sit down while showering.</p> <p>VALUE IV shower seat. Fixed height. The seat is foldable.<br/>         Materials: Diecast aluminium, Rilsan-coated stainless steel and vacuum-moulded polyethylene (PE).<br/>         Incl. fixing bolts for brick wall. Plugs depend on wall type and must be ordered separately.</p>                                                                                                                                                                                                                         |
| <b>Adjustability</b><br><b>Product capacity</b>        | Fixed mounting<br>Maximum user mass: 135 kg.<br>When the seat is loaded with 135 kg, each wall fastener is loaded with 125 kg.<br><i>The seat has been tested statically with safety factor 1.5, which is a requirement according to the applied assistive products standards.</i>                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Materials</b><br><b>Ergonomy</b><br><b>Features</b> | <p>The wall profile is made of diecast aluminium, and the frame and the leg is made of Rilsan-coated stainless steel. Rilsan provides good protection against corrosion and is comfortable to touch and easy to clean. The seat is made of vacuum-moulded polyethylene (PE).</p> <p>The anti-slip pattern of the seat ensures drainage of water in both the horizontal and the vertical position.<br/>         The seat is slightly curved towards the middle, which helps the user stay in the middle of the seat.</p> <p>The folding shower seat is designed for installation in the shower area and can be easily folded up and down as required.<br/>         The seat and backrest can be easily removed for cleaning.</p> |



|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b>     | VALUE IV folding shower seat with supporting leg, fixed height                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Item number</b> | R1622000112                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Warranty</b>    | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Cleaning</b>    | <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p><b>DISINFECTANTS:</b><br/>Clinimax Difficil-S<br/>Erisan OXY+<br/>Sprit (husholdningssprit)<br/>Suma Bac D10<br/>Virkon<br/>Rodalon Desinfektionsspray<br/>Lysoformin Spezial<br/>CaviCide<br/>Sufanios Premium<br/>Tvättsprit 70%<br/>WipeClean UniversalCleaning<br/>Chlorine Disinfection Wipes<br/>Klorilli (Kiilto Clean)</p> <p><b>NON-DISINFECTANTS:</b><br/>Suma Star free D1<br/>Taski 100<br/>Taski 200<br/>Milizid Citro</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p> |



# PRESSALIT®

|                    |                                                                |
|--------------------|----------------------------------------------------------------|
| <b>Product</b>     | VALUE IV folding shower seat with supporting leg, fixed height |
| <b>Item number</b> | R1622000112                                                    |

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®

13. juli 2026