



# PRESSALIT®

|                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b><br><b>Item number</b> | Changing table mattress 1800 mm<br>R8704000                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Colour</b>                        | VVS: 777630760<br>RSK: 8565134<br>EAN: 5708590366098<br>HMI: 125811<br><br>White<br>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Product description</b>           | Mattresses for daycare changing tables are available in several colors and in lengths that are adapted to the different changing table lengths.<br><br>Accessories for changing tables.<br>Mattress 760 mm x 1800 mm x 30 mm, for 1800 mm changing table.<br>Phthalate-free.                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Materials</b>                     | Cover: Techmadin, phthalate-free<br>Padding: Foam                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Cleaning</b>                      | <b>Upholstery (Cover):</b> Techmadin, phthalate-free<br>Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and wash regularly at 60-90°, without fabric softener. Shrinkage after washing may occur. Do not turn inside out before washing. Can be tumble-dried at max. 60°. Can be disinfected with antiseptic wipes. Alcohol % < 70.<br>Do not use disinfectants intended for hard surfaces.<br><br><b>Padding: (foam)</b><br>The foam contains only DINP softeners which are approved by the EU for the use in articles for children.<br>Do not wash the foam.<br>Do not use the mattress without the cover for reasons of hygiene. |



# PRESSALIT®

|                    |                                 |
|--------------------|---------------------------------|
| <b>Product</b>     | Changing table mattress 1800 mm |
| <b>Item number</b> | R8704000                        |

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®

9. september 2026