



PRESSALIT®

Product Item number	Pressalit Public Premium Toilet paper holder for 1 roll of toilet paper + 1 spare roll, polished steel PX10210
Product description	VVS: 776130240 NRF: 1438648 RSK: 8920358 EAN: 5708590392073 Toilet paper holder in steel housing with space for spare roll. For wall mounting. The holder is lockable and includes mounting kit.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com