



# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Premium Soap Dispenser, 200 ml, brushed steel PX30710
<b>Product description</b>	VVS:774684215 RSK:8920370 EAN:5708590392172  Stainless steel soap dispenser for wall mounting. One-hand operated and with drip-free soap pump with a 200 ml container. Designed for regular liquid soap or hand lotion. Incl. mounting kit.
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)