




# PRESSALIT®

<b>Product</b>	Mixer tap
<b>Item number</b>	RT633
<b>Product description</b>	<p>EAN: 5708590396309</p>  <p>Mixer tap for daycare changing tables and basin vanity unit</p>
<b>Materials</b>	Brass with chrome plating
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
<b>Cleaning</b>	The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

[pressalit.com](http://pressalit.com)